
3i T3 IMPLANT™

PRESERVATION
BY DESIGN™

BIOMET 3i™
PROVIDING SOLUTIONS - ONE PATIENT AT A TIME™

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
1. Lazzara RJ,[†] Porter SS. Platform Switching: A New Concept In Implant Dentistry For Controlling Post Restorative Crestal Bone Levels. *Int J Periodontics Restorative Dent.* 2006;26:9-17.
2. Vela X[†], Méndez V[†], Rodríguez X[†], Segalá M[†], Tarnow DP[†]. Crestal Bone Changes On Platform-Switched Implants And Adjacent Teeth When The Tooth Implant Distance Is Less Than 1.5mm. *Int J Periodontics Restorative Dent.* 2012 Apr;32(2):149-155.
3. Fransson C, Lekholm U, Jemt T, Berglundh T. Prevalence Of Subjects With Progressive Bone Loss At Implants. *Clinical Oral Implants Research.* 2005;16:440-446.
4. Zitzmann NU, Berglundh T. Definition And Prevalence Of Peri-Implant Diseases. *Journal of Clinical Periodontology.* 2008;35:286-291.
5. Sullivan DY, Sherwood RL, Porter SS. Long-Term Performance Of Osseotite® Implants: A 6-Year Follow-Up. *Compendium.* 2001 Apr;Vol.22, No.4.
6. Javed F, Romanos GE. The Role Of Primary Stability For Successful Immediate Loading Of Dental Implants. A Literature Review. *J Dent.* 2010 Aug;38(8):612-20. Epub 2010 Jun 11. Review.

[†]Dr Goené, Dr Lazzara, Dr Makigusa, Dr Méndez, Dr Rodríguez, Dr Segalá, Dr Tarnow and Dr Vela have financial relationships with BIOMET 3i LLC resulting from speaking engagements, consulting engagements and other retained services.

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
7. Testori T†. Clinical Factors Related To Implant Stability With Tapered Implants. *J Implant And Res Dent*. Vol. 2, No. 1, 2010.

8. Östman PO†, Wennerberg A, Ekström A, et al. Immediate Occlusal Loading Of NanoTite™ Tapered Implants: A Prospective 1-Year Clinical And Radiographic Study. *Clin Implant Dent Relat Res* 2012 Jan 17. [Epub ahead of print]

9. Block MS†. Placement Of Implants Into Fresh Molar Sites: Results Of 35 Cases. *J Oral Maxillofac Surg*. 2011 Jan;69(1):170-4.

10. Lin A, Wang CJ, Kelly J, Gubbi P†, Nishimura I. The Role Of Titanium Implant Surface Modification With Hydroxyapatite Nanoparticles In Progressive Early Bone-Implant Fixation In Vivo. *Int J Oral Maxillofac Implants* 2009 Sep-Oct;24(5):808-816.

11. Davies JE. Bone bonding at natural and biomaterial surfaces. *Biomaterials*. 2007 Dec;28(34):5058-67.

12. Mendes VC, Moineddin R, Davies JE. The effect of discrete calcium phosphate nanocrystals on bone bonding to titanium surfaces. *Biomaterials*. 2007;28:4748-4755.

13. Kikuchi L, Park JY, Victor C, Davies JE. Platelet interactions with calcium-phosphate-coated surfaces. *Biomaterials* 2005 Sep;26(26):5285-95.

14. Park JY, Gemmell CH, Davies JE. Platelet interactions with titanium: modulation of platelet activity by surface topography. *Biomaterials* 2001 Oct;22(19):2671-82.

15. Kuzysk PR, Schemitsch EH. The basic science of peri-implant bone healing. *Indian J Orthop*. 2011 Mar;45(2):108-15.

16. Davies JE. Understanding Peri-Implant Endosseous Healing. *J Dent Educ*. 2003 Aug;67(8):932-49

†Dr Block, Dr Gubbi, Dr Östman and Dr Testori have financial relationships with BIOMET 3i LLC resulting from speaking engagements, consulting engagements and other retained services.

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

Zetterqvist et al. A Prospective, Multicenter, Randomized Controlled 5-Year Study Of Hybrid And Fully Etched Implants For The Incidence Of Peri-implantitis. *J Periodontol* April 2010.

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17. Zetterqvist et al. A Prospective, Multicenter, Randomized Controlled 5-Year Study Of Hybrid And Fully Etched Implants For The Incidence Of Peri-implantitis. *J Periodontol* April 2010.
18. Wennerberg A, Albrektsson T. On implant surfaces: a review of current knowledge and opinions. *Int J Oral Maxillofac Implants*. 2010 Jan-Feb;25(1):63-74. Review.
19. Gubbi P, Towse R. Quantitative and Qualitative Characterization of Various Dental Implant Surfaces. Presented at the European Association of Osseointegration, 20th Annual Scientific Meeting; October 2012; Copenhagen, Denmark.

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

20. Baumgarten H¹, Meltzer A¹. Improving Outcomes While Employing Accelerated Treatment Protocols Within The Aesthetic Zone: From Single Tooth To Full Arch Restorations. Presented at the Academy of Osseointegration, 27th Annual Meeting; March 2012; Phoenix, AZ.

21. Lazzara R.¹ Dental Implant System Design and the Potential Impact on Long-Term Aesthetics: A Review of the 3i T3™ Tapered Implant, ART1193EU 3i T3 White Paper. BIOMET 3i, Palm Beach Gardens, Florida, USA.

22. Suttin Z, Towse R, Cruz J. A Novel Method for Assessing Implant-Abutment Connection Seal Robustness. BIOMET 3i, Palm Beach Gardens, Florida, USA. Poster Presentation, Academy of Osseointegration, 27th Annual Meeting; March 2012; Phoenix, AZ. http://biomet3i.com/Pdf/Posters/Poster_Seal%20Study_ZS_AO2012_no%20logo.pdf. Testing done by BIOMET 3i, Palm Beach Gardens, FL; n = 20.

23. Al-Jadaa A, Attin T, Peltomaki T, Schmidlin PR Comparison of three in vitro implant leakage testing methods. Clinica Oral Implant Research 2013 Dec. 16. doi: 10.1111/clr/12314

24. Byrne D, Jacobs S, O'Connell B, Houston F, Claffey N. Preloads generated with repeated tightening in three types of screws used in dental implant assemblies. J. Prosthodont. 2006 May-Jun;15(3):164-71.

¹Dr Baumgarten, Dr Lazzara, Dr Meltzer and Dr Östman have financial relationships with BIOMET 3i LLC resulting from speaking engagements, consulting engagements and other retained services.

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²⁵ Boitel N, Andreoni C, Grunder U[†], Naef R, Meyenberg K[†]. A Three Year Prospective, Multicenter, Randomized-Controlled Study Evaluating Platform-Switching For The Preservation Of Peri-Implant Bone Levels. *Poster Presentation P83: Academy of Osseointegration, 26th Annual Meeting*; 2011 March 3-5; Washington DC.

²⁶ Östman PO,[†] Wennerberg A, Albrektsson T. Immediate Occlusal Loading Of Nanotite Preval Implants: A Prospective 1-Year Clinical And Radiographic Study. *Clin Implant Dent Relat Res*. 2010 Mar;12(1):39-47.

[†]Dr Grunder, Dr Östman, Dr Meyenberg, Dr Scutella, Prof. Testori and Dr Vela have financial relationships with BIOMET 3i LLC resulting from speaking engagements, consulting engagements and other retained services.

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Encode, GingiHue, Gold-Tite, OSSEOTITE, PreFormance and PREVAIL are registered trademarks of BIOMET 3i LLC. BellaTek design, Providing Solutions - One Patient At A Time, 3i T3 and 3i T3 Implant design are trademarks of BIOMET 3i LLC. Endobon is a registered trademark of BIOMET Deutschland. ©2012 BIOMET 3i LLC.